


Equagel Cushions

General Cushion

- For moderate cushioning needs and prevention of pressure sores
- The thickness of the cushion is 2 inches
- User Weight 150kg



Protector Cushion

- For advanced cushioning needs
 The thickness of the cushion is 2.5 inches
 - User Weight 147kg



Adjustable Protector Cushion

- For customisable cushioning needs
- The thickness of the cushion is 2.5 inches
- Two different pelvic contour inserts are included to increase the therapeutic gel under the pelvis
- User Weight 160kg

